

Pearlescent NanoWhite (35%) In-Office Kit

Professional Directions for Use

Overview

Pearlescent NanoWhite is a professional 35% in-office teeth-whitening system. Its nanotechnology gel penetrates the enamel to target stubborn stains and discoloration, and is suitable for both vital and non-vital teeth whitening.

Product details

Concentration: 35% (in-office)

Application: In-office / chairside, dental professional supervised

Treatment time: From 18 minutes per cycle (1 cycle)

Suitable for: Vital and non-vital teeth

Technology: Nanotechnology whitening gel

Key points

- Effective, visible results in as little as one 18-minute session.
- Specially formulated concentration helps minimise sensitivity while protecting enamel.
- Administered by a trusted dental professional in a clinic setting.
- Simple post-whitening care routine helps maintain results.

Innobrace Orthodontics - innobrace.com