

Pearlescent nanoWhite Teeth Whitening

A New Standard in Professional Teeth Whitening

Pearlescent Nanowhite represents a breakthrough in in-office whitening. By combining the power of titanium nanoparticles with advanced violet light technology, it delivers faster, more effective results while using a lower concentration of peroxides—enhancing both safety and patient comfort. Patients can achieve visibly whiter smiles in less time, with a treatment experience designed to be gentle, efficient, and reliable.

Why Choose Pearlescent Nanowhite?

Guaranteed Safety

Formulated with a lower concentration of peroxides to help preserve enamel while delivering effective whitening results.

Zero Sensitivity

Designed for maximum comfort before, during, and after treatment, ensuring a pleasant patient experience.

Superior Performance

Achieves noticeable results in less time, helping optimise clinical efficiency and reduce chair time. When used under the care of a dental professional, Pearlescent NanoWhite offers effective results that can usually be seen in just one session.


Quick & Fuss-free

In office treatment can be completed in less than an hour for both upper and lower arches.



Pearlescent
nanoWhite
Because every smile deserves to
shine with Pearlescent NanoWhite!

Authorised Dental Professional



 **Innobrace
Orthodontics**
we give the world a smile

Innobrace Orthodontics Pte Ltd

28 Senang Crescent, #02-01/02,
Bizhub 28, Singapore 416601.

 **(65) 6446 6317**
 marketing@innobrace.com
 www.innobrace.com

Call / Whatsapp :
 **(65) 9117 1158, (65) 9176 6746**


 **Innobrace
Orthodontics**
we give the world a smile

Pearlescent
nanoWhite

25% Hydrogen
Peroxide

Dental Professional Kit

Unlock Your
Radiant Smile!



Pearlescent nanoWhite (25%) Dental Professional Kit



What causes tooth discoloration?

Some common causes of tooth discoloration include ageing, smoking, ingestion of staining substances over a period of time, eg. Coffee, tea, red wine, curries, tobacco, colas may also contribute to tooth discoloration.

Is tooth whitening safe?

Tooth whitening is an effective, traditional, safe and easy technique that will please patients with its final natural-looking results. Over the years, manufacturers have greatly improved and advanced the formulation needed to achieve desired results while ensuring patients' comfort and safety. It can be considered as a first treatment option for cases of color alteration because it is a low cost technique with little side effects and does not imply wear of dental structure.

Am I suitable for teeth whitening?

Almost anyone. However, as results may differ between individuals based on your medical and dental history, dental professionals will be able to assess every patient's teeth condition and advise if you are suitable to undergo whitening treatment or they may tailor the treatment plan to fit individual needs.

FREQUENTLY ASKED QUESTIONS

Q: Why Choose Pearlescent NanoWhite?

A: Pearlescent's innovative Nanotechnology formula is a whitening gel that penetrates deep into the enamel, effectively targeting stubborn stains and discoloration for a brighter smile. Pearlescent Nano-white can be used for vital and non-vital teeth whitening.

Q: How long is a whitening session?

A: Most patients achieve effective **visible results in as little as 18mins treatment time (1 cycle)**, leaving you with a confident, radiant smile.

Q: Is there any discomfort during the treatment?

A: Gentle on your teeth: Experience little to no sensitivity as we protect teeth enamel from harsh whitening concentration yet deliver effective visible results in as little as **1 session**.

Safe & Comfortable: Pearlescent Nanowhite is specially formulated to minimize sensitivity, our kit ensures a gentle yet effective whitening experience for both vital and non-vital teeth. Be assured that you are in safe hands of your trusted dental professionals as they guide you through your entire whitening experience in the comfort of a clinic setting.

Q: Is whitening treatment expensive?

A: One of the Best Cost benefit in the market: Pearlescent offers one of the lowest treatment costs in the market with effective visible whitening results.

Q: How long do whitening results last?

A: Keeping to simple post whitening care routine helps in the maintenance of your whitening treatment, which differ from individuals to individuals.

Indications for Use

Pearlescent Nanowhite is ideal for advanced teeth whitening in:

- Patients seeking quick and effective results
- Professionals aiming to streamline procedures and enhance patient satisfaction

How It Works

1. Consultation:

Your dental professional will assess your smile and recommend a personalized treatment plan. Your teeth will be prepped and cleaned prior to a whitening session.

2. Application:

The whitening gel is applied to your teeth directly. Special care will be taken by your doctor to ensure a comfortable experience.

3. Results:

Enjoy your brighter smile in as little as 1 treatment cycle (18minutes)! Follow your dental professional's post-treatment care tips/advice to maintain your results.

Perfect for Everyone

Whether you're preparing for a special occasion or just want to enhance your everyday smile, the Pearlescent NanoWhite Kit is a perfect choice for your budget and lifestyle!

Ready to Shine?

Ask your dental professional about the **Pearlescent NanoWhite Teeth Whitening Kit** today and take the first step toward a dazzling smile!

